

# SPACEPLUS & ONLINESUITEPLUS SPACEPLUS & ONLINESUITEPLUS

ADVANCED - SUPERIOR - DIGITAL



# We protect and improve the health of people around the world.

**OUR VISION** 



#### B. BRAUN AIS – 70 YEARS OF EXPERTISE





#### B. BRAUN – 70 YEARS OF EXPERTISE



Infusomat' fm



#### B. BRAUN – 70 YEARS OF EXPERTISE



# Trust in experience Quality of life



# **ADVANCED** 0-+-0---+00



6













**INTUITIVE USER INTERFACE** 



AT A GLANCE



DOSEGUARD -

THE SAFER, THE

BETTER

DRUG COLOR

CODING(ISO 26825:2020)

-11			



VARIOUS THERAPY MODES













# **OnlineSuite**plus







...central digital platform

...part of smart infusion therapy

# The next generation of integrated infusion pump systems





Not just an infusion pump!

# OnlineSuite<sup>plus</sup> App Store



. . . . . . . .

# Platform Administration





# **DeviceManager**<sup>plus</sup>



# Make Progress Not Miles



Smart Asset Management Track your device and know all details!



Swift & Simple Data Management Software, Drug library, configurations are always up to date!



#### B. Braun Health Cloud

Application is fully integrated into the B. Braun Health Cloud for instant access to new updates!



# **B. Braun Health Cloud Connector**



# Smart Technologies

Smart Hospitals



# BBRAUN



#### Central system

Integrate B. Braun medical devices into hospitals processes and workflows !



#### Access to latest updates!

With B. Braun Health Cloud you always access to the latest software updates for new functionalities or cybersecurity.



#### Fully digital process!

- Mail notification of the new update
- Automatically availabe in OnlineSuite<sup>plus</sup>
- Online distribution via DeviceManager<sup>plus</sup>





# DoseLinkplus



# Smart Technologies Smart Hospitals



HealthCare

ascom

Allscripts Gerner



**Central application for system interoperabilty** DoseLink<sup>plus</sup> provides integration with other systems in a standarized HL7 interface, acc. to IHE standards



#### Secure interoperability

Pumps communicate with the platform with full encryption of data and acc. to digital certificates.



capsule

a Qualcomm Life company

Epic

GE Healthcare

Dräger

#### Seamless documentation

Therapy data is seamlessly documented in the EHR and the pumps are equipped with data buffering

Drug Library Manager<sup>plus</sup>



# **Even more performance** Even more Space



## **BBRAUN**



## DrugLibrary Manager

Administrate and configure your drug library centrally!



#### The safer. The better.

Reduces the risk of medication errors. Various drugs (max 10.000), drug categories (up to 31), patient profiles (up to 32) and specific care units (up to 31) can be programmed.



#### Drug color labels

International standard colors are available for different drug groups and can be flexibly assigned.



#### Therapy Modes – keep control.

All therapy modes are unlocked for the pump and visualized on the screen and can thus be monitored and controlled in a more targeted manner.





# DoseTracplus



# No longer searching

## any report data



#### **Display history infusion data** Wide range of data to help avoid pot

Wide range of data to help avoid potential medical errors!



#### **Reporting & Analytics**

Retrospective reporting reaching the highest level of patient safety to prevent medication errors.

- Drug library compliance
- Alarm frequency
- Soft & hard limit notifications
- Bolus reports



#### Drug database optimization

Creation of the drug database and ensuring correct application based on clinical practice.



## **BBRAUN**





# **OneView**<sup>plus</sup>



# Interact with technology Interface with life





## Smart Alarm Management

Display on monitors and on mobile phones. Centrally informed about every alarms.



#### **Silent Hospital**

Keeps staff informed without audible alarms in patient room. Less stress!





E NeurolCU 🔘		=	NeuroIC	J
Filter: Night Shift			Filter: Night S	hift
Room 201				21
201_1		Drug Insulin 1. Bed Bed 301_3	.0 mg/ml <sup>Serial Number</sup> 87523	
Room 301				161
301_1 (()) (()) (()) (()) (()) (()) (()) (	$\Leftrightarrow$	Adrenalii Bed Bed 201_1	n4 0.001 mE <sup>Serial Number</sup> 349258	
Room 302		201_1		
		1 1		293 1
302_1 (()) 302_2 ×		<sup>вед</sup> Bed 301_1	Serial Number 111197	
 ⊘6 ∧2 ⊗4		<b>⊳</b> 6	<u>^</u>	
Rooms Alarms Remaining Times		Rooms	Alarms	Re
		_	_	-

2 Minutes

16 Minutes

293 Minutes

Pump type

Pump type

Pump type

N)5

Remaining Times

# The next generation of integrated infusion pump systems



Supported by an innovative infusion pump system







# **Compact and flexible concept** Space<sup>plus</sup> weights only 1,9 kg and has an integrated handle for flexbile use of the system.



#### **Resistive Touchscreen & Handling** Get in touch with a new and intuitive therapy experience. Operating feels as familiar as using your smartphone. Even with wet gloves.



#### Therapy Modes – keep control.

All therapy modes are visualized on the screen and can thus be monitored and controlled in a more targeted manner.





# Spaceplus



#### Same UI interface

To reduce training effort, Space<sup>plus</sup> has the same user interface logic for volumetric and syringe pump.



#### No physical connectors

Spaceplus communicates with the Station via infrared. Corrosion of physical connectors are completely avoided.





#### Long lasting battery

Long lasting batteries and an integrated power supply liberate care and set the staff free from inflexible wall sockets.



## Station Space<sup>plus</sup>



#### Superior workplace design

The benchmark in ergonomy: The modular system enables configurations of up to 24 pumps. Its portability mobilizes therapy.



#### Central audible & visual alarm

The Space<sup>plus</sup> Cover centrally provides the audible and visual alarm of the pumps. Integrated brightness sensor automatically adjusts acc. to the environment.

#### **Optimized fixation**

The Station clamp is centrally placed to better bring the center of gravity closer to the IV pole. Standard mounting for vertical installation, optionally for horizontal mounting for DIN rails available.









## Connectivity Spaceplus



#### Integrated WiFi

Experience cordless care. Pumps can be connected simply via integrated WiFi or – when stacked in a Station – via Ethernet.



### LAN – Ethernet connection Optionally, a data module can be connected

at the bottom of the Station to set up a LAN connection.



#### **Communication to the Data module**

Space<sup>plus</sup> pumps automatically recognize the data module to set up the communication automatically.



#### HDMI

Data module is equipped with HDMI connection for future-oriented connection options.







## Design your workplace



#### Stand-alone - Transport

Space<sup>plus</sup> pumps can be easily stacked to each other for attachment to an IV pole or for transport.





#### **Docking station**

Space<sup>plus</sup> pumps can be inserted and taken out of the Station w/o interruption of running therapy.







# The Space<sup>plus</sup> MRI Station supports an infusion therapy for one patient at a time during MRI examinations w/o artifacts.

MRI

## Accurate & Safe



#### **IP44 Moisture Protection**

The pump is protected from moisture and splashing water – no matter what

#### **Delivery accuracy**

- Syringe pump +/- 2%
- Volumetric pump +/- 3%



#### Freeflow is realiably prevented

The unique and patented syringe piston brake prevents freeflow and handling mistakes during syringe change.





# Space<sup>plus</sup> Perfusor – Inserting syringe





## Disposable Portfolio



#### In-house production

Infusion therapy has been one of the focus areas of B.Braun for > 75 yrs and has all production sites worldwide spread under own high quality management requirements.



A wide range of dedicated disposables Standard, SafeSet, Light protected, NRFit, ENFit, Flushing Set, Transfusion 2 mL – 50/60 mL syringes



#### Sustainable & safe raw materials Safe raw materials are always used for the production of our systems, in which

sustainability is also one of our priorities.







Enterale infusion set

Standard IV-set

Oncology infusion set



## Disposable Portfolio





## Disposable Portfolio




Disposable Portfolio





## Disposable Portfolio





## Disposable Portfolio



# Supports all departments within your hospital



















#### Superior workplace design

The benchmark in ergonomy: The modular system enables configurations of up to 24 pumps. Its portability mobilizes therapy.









#### Superior workplace design

The benchmark in ergonomic: The modular system enables configurations of up to 24 pumps. Its portability mobilizes therapy.



#### **Touchscreen & Handling**

Intuitive user interface with all essential therapy data at a glance, including drug colors labels.











#### Superior workplace design

The benchmark in ergonomy: The modular system enables configurations of up to 24 pumps. Its portability mobilizes therapy.



#### Touchscreen & Handling

Intuitive user interface with all essential therapy data at a glance, including drug colors labels.



#### DoseGuard – The safer. The better.

DoseGuard reduces the risk of medication errors. Various drugs, drug categories and specific care units can be programmed.







## Alarm Management



#### Drug name always visible

During therapy, in standby or in case of alarm? Yes, the drug name is visible at all times.



#### Smart Alarm Management

Every alarm with the cause is displayed on monitors or on mobile phones.



#### Silent hospital approach

Space<sup>plus</sup> makes it possible. Alarms can be suppressed in case of confirmed distributed alarm management.







# Alarm confirmation

\_\_\_\_

≡	NeurolCU 🥘		
Filter: Night Shift			
Roc	om 201		
20	1_1		
Room 301			
1000	1_1 301_2 301_3		
Room 302			
	2_1 302_2 ×		
R	ooms Alarms Remaining Times		



## Accurate Medication Management



#### Precise drive & peristalsis

Rethink rigor in drug administration. Both, Perfusor and Infusomat feature a totally redesigned drive and peristalsis technology.



#### Administration of Inotrops – Take Over Mode

Reduce interuption of delivery of drugs with short half-life with the automatic syringe change to keep the blood pressure constant.



#### Seamless therapy data into EMR

Connection to EMR are seamless via standardized HL7 / IHE interface. Wifi and data buffering makes the system very flexible.









**One pump for all therapies** Continuous, Enteral, TCI or PCA? Space<sup>plus</sup> supports all therapies.



#### **Reduce workload**

Synchronization of patient related data within the Station.



#### Seamless therapy data into EMR

Connection to EMR are seamless via standardized HL7 / IHE interface. Wifi and data buffering makes the system very flexible.





## Target Controlled Infusion



#### Multiple algorithms available

Space<sup>plus</sup> supports a wide range of algorithms that are fully configurable in the Drug LibraryManager<sup>plus</sup> to reduce handling steps.



#### **Graphical visualization**

The large color screen displays the therapy in full graphics. Titration of the target can be done very quickly and is immediately graphically visible.



Awake – one click to wake up the patient Terminate anaesthesia? Press "awake" and the target will automatically be set to 0. Concentration reduction remains graphically visibile with wake-up time indication.





## Patient Controlled Analgesia



#### **PCA** profiles

Space<sup>plus</sup> supports PCA, PCEA, IPCA and PIEB. All profiles are to configured in the Drug LibraryManager<sup>plus</sup>.



#### **Graphical visualization**

The large color screen displays the A/D ratio in full graphical form.



#### Lockbox

Space<sup>plus</sup> has a PCA lockbox to prevent tampering the drugs.



#### **Data lock**

Extra data locks to give respective access / rights to special pain team.











#### **Small & Intuitive**

Space<sup>plus</sup> is very small and can be easily attached to a vertical / horizontal IV pole / DIN rail.



#### **Medication safety**

DoseGuard<sup>™</sup> is the function to define limit ranges to reduce medication errors. Next step is the possibility of closed loop medication management.



#### Silent hospital approach

To reduce alarm fatique as well as to give patient the necessary rest during administration of chemotherapy.





Smart Chemotherapy Management

# PREVENT MEDICATION ERRORS



# REASONS FOR MEDICATION ERRORS



88% of the errors were detected in the day hospital while 12% reached the patient.<sup>3</sup>

1 Ford et al (2006): Study of Medication Errors on a Community Hospital Oncology Ward. Journal of Oncology Practice, 2006, 2 (4), 149-154. available at: https://ascopubs.org/doi/full/10.1200/jop.2006.2.4.149; accessed 05-06-2019.

2 Schulmeister L (1999): Chemotherapy medication errors: descriptions, severity, and contributing factors. Oncol Nurs Forum. 1999 26(6): 1033-42

3 Villar J, et al. Errores asociados con la prescripción, validación, preparación y administración de medicamentos citostáticos. Farmacia hospitalaria: órgano oficial de expresión científica de la Sociedad Española de Farmacia Hospitalaria. 2008; 32(3):163-169.



## Smart Chemotherapy Management



**Organizational Measures** 



## Smart Chemotherapy Management



#### **Closed Loop Medication Management** Space<sup>plus</sup> supports bi-directional communication (AutoProgramming), meaning

not only sends out data for documentation but also receives therapy order.



#### **AutoProgramming**

Digital process to support barcode scanning for therapy validation, patient identification, therapy related infusion data and documentation.



#### Partners in smart chemotherapy

Already with the Space generation, there were close collaboration with various partners.







## Medication Management



#### **DrugLibrary Manager**

Administrate your medication data centrally and configure your drug library adapted to every ward / even patient profile. Collaboration platform for hospital medication safety team



#### Configuration settings per care area

Certain departments, medication or therapies require specific settings regarding pressure, KVO rates, patient limits etc. This can be managed centrally.







## Medication Management



#### **DrugLibrary Manager**

Administrate your medication data centrally and configure your drug library adapted to every ward / even patient profile. Collaboration platform for hospital medication safety team



#### Configuration settings per care area

Certain departments, medication or therapies require specific settings regarding pressure, KVO rates, patient limits etc. This can be managed centrally.



#### Automatic validation process

Adding or adjusting a medication setting is fully supported by a validation process. Every change is documented.

2 3 5 **Care Area Drug Data** Concentrations Therapies **Care Areas** Variant Modify therapy to **Define therapy** Define general Add multiple Assign drug to care area specific setting and Drug settings Care Areas concentration patient profiles requirements

> Unlimited change / release history per drug





Automatic notifications for released / ready DrugLists



Change



## Medication Management



#### **DrugLibrary Manager**

Administrate your medication data centrally and configure your drug library adapted to every ward / even patient profile. Collaboration platform for hospital medication safety team



#### Configuration settings per care area

Certain departments, medication or therapies require specific settings regarding pressure, KVO rates, patient limits etc. This can be managed centrally.



<u>∎</u> ©©0

**B**BRAUN

#### Automatic validation process

Adding or adjusting a medication setting is fully supported by a validation process. Every change is documented.

#### **Compliance reporting**

DoseTrac<sup>plus</sup> provides insights into how the drug library is used in order to optimize the drug library.







# BIOMEDICAL DEPARTMENT

TAT

Ш

Add drug

WA HAT IS NOT STATE

Contraction of the second seco

Terrap Step

Crute drug Roury

0



## System construction



#### **Reliable and rigid Construction** Oops! Fallen? The pump – especially the drive arm - is designed to withstand the challenges of everyday life.







- Electronics stored underneath top cover
- Integrated Wi-Fi
- No soldered connections to accesory plug
- New drive head & drive housing (4)

- Integrated power supply integrated
- New and optimized aluminium drive mechanism
- Improved sealing

3

Improved contour & sealing of drive arm (4) 61

B. Braun Melsungen AG



## System construction



#### Reliable and rigid Construction

Oops! Fallen? The pump – especially the drive arm – is designed to withstand the challenges of everyday life.



#### **IP44 Moisture Protection**

The pump is protected from moisture and splashing water – no matter what.





## System construction



#### Reliable and rigid Construction

Oops! Fallen? The pump – especially the drive arm – is designed to withstand the challenges of everyday life.



### IP44 Moisture Protection

The pump is protected from moisture and splashing water – no matter what.



#### Efficient technical safety check

TSC is required every 24 months w/o opening the pump and preventive replacement of parts (Battery > 6 yrs). User will automatically get a pop-up.





TSC information on pump Pop-up information





## Make progress – No miles



#### **Central management**

With DevicerManager<sup>plus</sup> you will manage all infusion pumps in a central digital platform.



Swift & Simple Data Management Software, Drug library, configurations are always up to date!



#### Track & Trace

Localization is an important part where the DeviceManager<sup>plus</sup> will support you with Wifi triangulation.







## System Architecture

OnlineSuite<sup>plus</sup> Apps



Cybersecurity







#### **Data encryption**

Communication between Space<sup>plus</sup> and OnlineSuite<sup>plus</sup> is fully encrypted and secured with digital signatures



## Swift & Simple Data Management

Software is always up to date! Drug Library are continuously updated according to hospital guidelines w/o interrupting clinical workflows.



**24/7 access to latest updates & upgrades** With the B.Braun Health Cloud you will receive notifications of new available updates and they are immediately available. Direct Access to Health Cloud Software updates, upgrades & Cybersecurity fixes



**OnlineSuite**<sup>plus</sup>

signatures





## System Requirements



#### **Operating System**

- Linux (RHEL, Debian, Ubuntu,...)
- Windows Server 2019

•••	
•••	
•••	

#### Server Hardware

- CPU: 4x 2.7 Ghz
- RAM: 16 GB
- HDD: 100 GB







#### Installation Opportunities В С Α OnlineSuite<sup>plus</sup> OnlineSuite<sup>plus</sup> OnlineSuite<sup>plus</sup> **B BRAUN B BRAUN B BRAUN** Docker Runtime CE\* Docker Runtime CE\* Docker Runtime CE\* Linux Linux Linux (Debian, Ubuntu, Windows Server 2019 / Host OS with Hypervisor RedHat EL etc.) Windows Server 2022 (VMware ESX, etc.) (Hyper-V) Server Hardware Server Hardware Server Hardware

